


















SERVICE SOLUTIONS
**Compressor
 Rebuild Kits**



APPLICATION	COMPRESSOR	DESCRIPTION	OXUS PART #
Invacare Perfecto2V DeVilbiss 5 Liter	 GSE ZW280	 GSE 280 Top End Rebuild Kit	RS-01001 RS-01018 (220V - Compatible in Europe)
Invacare Platinum 5 Invacare PerfectO2 Airsep New Life Nidek Nuvo 8 and 10	 THOMAS 2650/2660	 THOMAS 2650/2660 Rebuild Kit	RS-01003 (Top End Rebuild Kit) RS-01004 (Top End Rebuild Kit w/ Cylinders)
Respironics Everflo Nidek Nuvo Lite	 THOMAS 2450	 THOMAS 2450 Top End Rebuild Kit w/ Cylinders	RS-01005
Invacare Platinum 5	 GAST 82R	 Gast 82R Top End Rebuild Kit w/ Cylinders	RS-01006
Puritan Bennett	 GAST 72R	 Gast 72R Top End Rebuild Kit w/ Cylinders	RS-01007
Invacare Platinum 10 (OEM Kit # 1188159)	 GSE ZW400-D2	 GSE ZW400-D2 Top End Rebuild Kit w/ Cylinders	RS-01015
Invacare Platinum 10 (OEM Kit # 1153687)	 THOMAS 2660	 THOMAS 2660 Top End Rebuild Kit w/ Cylinders (≥ 0.50" stroke)	RS-01017

Oxus America
 2046 Brown Rd.
 Auburn Hills, MI 48326
 Phone: 248-475-0925 | Toll Free: 888-475-1568
 sales@oxusamerica.com
 www.oxusamerica.com



Oxus is an engineering, manufacturing, and service company specializing in medical devices and gas separation technology. We are a FDA registered, ISO 13485 Certified OEM Manufacturer and have been in business for over 20 years.