



SERVICE SOLUTIONS

POC Sieve Replacements



Oxus has developed a high quality aftermarket solution for POC replacement sieves with an eye toward cost reduction. Oxus sieves are **100% functionally checked** prior to shipment and come with a **1 year warranty** against leaks and rupture.

Oxus replacement POC sieves include:

- Respironics SimplyGo
- AirSep Freestyle 3
- Inogen G3/OxyGo (4 setting and 5 setting models)
- Inogen G5/Rove 6/OxyGo Next

All Oxus POC replacement sieves:

- Engineered with high quality materials and construction to meet or exceed OEM specifications
- Use only factory new zeolite material
- Critical components are ultrasonically cleaned
- Made in the USA

Oxus America
2046 Brown Rd.
Auburn Hills, MI 48326
Phone: 248-475-0925 | Toll Free: 888-475-1568
sales@oxusamerica.com
www.oxusamerica.com



Oxus is an engineering, manufacturing, and service company specializing in medical devices and gas separation technology. We are a FDA registered, ISO 13485 Certified OEM Manufacturer and have been in business for over 20 years.

APPLICATION	OEM PART #	DESCRIPTION	OXUS PART #
 Respironics SimplyGo	1093443	 SimplyGo Sieve Asm.	RS-01010-1 (individual set) RS-01010-12 (12 set pack)
 AirSep Freestyle 3	BE227-1	 Freestyle 3 Sieve Asm.	RS-01011-1 (individual set)
 Inogen G3/OxyGo (Settings 1-4)	RP-320	 Inogen G3/OxyGo Sieve Asm. (low flow)	RS-01014-12 (12 set pack)
 Inogen G3/OxyGo (Settings 1-5)	RP-321	 Inogen G3/OxyGo Sieve Asm. (high flow)	RS-01016-12 (12 set pack)
 Inogen G5/Rove 6 OxyGo Next	RP-502	 Inogen G5/Rove 6/OxyGo Next Sieve Asm.	RS-01036-15 (15 pack)

Oxus America
 2046 Brown Rd.
 Auburn Hills, MI 48326
 Phone: 248-475-0925 | Toll Free: 888-475-1568
 sales@oxusamerica.com
 www.oxusamerica.com



Oxus is an engineering, manufacturing, and service company specializing in medical devices and gas separation technology. We are a FDA registered, ISO 13485 Certified OEM Manufacturer and have been in business for over 20 years.