



SERVICE SOLUTIONS

Stationary Concentrator Valves



All Oxus valves:

- Designed to meet or exceed OEM specifications
- 1 year warranty against wear and workmanship

Oxus 4-Way valves for Invacare:

- Completely refurbished including:
 - Complete disassembly and ultrasonic cleaning
 - Replacement of worn parts (including end caps and all seals)
 - Lubrication using a specific high performance and air-path compatible lubricant
- Fully tested to confirm valve performance prior to shipment
- Cores do not need to be sent prior to order fulfillment

Oxus Valves for Respiration EverFlo:

- 100 % new products
- Designed to eliminate body cracking

Oxus America

175 Kay Industrial Dr.

Orion Twp., MI 48359







Phone: 248-475-0925 | Toll Free: 888-475-1568

sales@oxusamerica.com

www.oxusamerica.com



Oxus is an engineering, manufacturing, and service company specializing in medical devices and gas separation technology. We are a FDA registered, MDSAP Certified, ISO 13485 Certified OEM Manufacturer and have been in business for over 20 years.

APPLICATION	OEM P/N	DESCRIPTION	OXUS P/N
All Invacare Concentrators	1101141	 4-Way Valve	RS-04000-1 (individual) RS-04000-24 (24 Pack)
Respironics EverFlo (Pre-June 2013)	1111248	 Valve Body - Internal Pilot - SMC Replacement	RS-04011-1 (individual) RS-04011-25 (25 Pack)
Respironics EverFlo (Pre-June 2013)	1114952	 Valve Body - Internal Pilot - ASCO Replacement	RS-04012-1 (individual) RS-04012-25 (25 Pack)
Respironics EverFlo (Post-June 2013)	1115003	 ASCO Valve Body Replacement (External Pilot)	RS-04013-1 (individual) RS-04013-25 (25 Pack)
All Respironics EverFlo Concentrators	1111249	 SMC Pilot Valve Asm.	RS-04004-1 (individual) RS-04004-25 (25 Pack)
All Respironics EverFlo (ASCO Applications)	1114954	 Pilot Valve Asm. (ASCO Replacement)	RS-04005-1 (individual) RS-04005-25 (25 Pack)

Oxus America
 175 Kay Industrial Dr.
 Orion Twp., MI 48359
 Phone: 248-475-0925 | Toll Free: 888-475-1568
sales@oxusamerica.com
www.oxusamerica.com



Oxus is an engineering, manufacturing, and service company specializing in medical devices and gas separation technology. We are a FDA registered, MDSAP Certified, ISO 13485 Certified OEM Manufacturer and have been in business for over 20 years.