



SERVICE SOLUTIONS

Thomas 2450 Compressors




Oxus can meet all of your Thomas 2450 Compressor needs with one of our multiple solutions to help increase your Everflo Concentrator maintenance capacity throughput:



Option 1: Oxus Thomas 2450 Refurbished Compressors

Oxus performs a complete Top End and Bottom End rebuild and each compressor is fully tested to meet or exceed OEM specifications, ensuring full compliance with FDA guidelines. Our refurbished compressors are available in 115V and 230V versions, can be shipped in bulk packs of 90 for the lowest shipping cost and come with a 1 year warranty on the Top End and 2 year warranty on the Bottom End.



Option 2: New Thomas 2450 Compressors

Don't have a core to send in? Purchasing a brand new Thomas 2450 Compressor from Oxus is a great solution when you do not have a core; no core return required. Our new compressors come with a 1 year warranty on the Top End and 2 year warranty on the Bottom End.

Thomas 2450 Compressor Accessory Kits


Increase efficiency and reduce bench time with Oxus' Thomas 2450 Compressor Accessory Kits.

Oxus America
175 Kay Industrial Dr.
Orion Twp., MI 48359
Phone: 248-475-0925 | Toll Free: 888-475-1568
sales@oxusamerica.com
www.oxusamerica.com



Oxus is an engineering, manufacturing, and service company specializing in medical devices and gas separation technology. We are a FDA registered, MDSAP Certified, ISO 13485 Certified OEM Manufacturer and have been in business for over 20 years.

Oxus Thomas 2450 Refurbished Compressor



1. Order by the Oxus Part # (order in packs of 90 for the lowest shipping cost)
115V Version - RS-03000-1 (Individual Pack) or RS-03000-90 (90 Pack)
230V Version - RS-03001-1 (Individual Pack) or RS-03001-90 (90 Pack)
2. Send Oxus your cores (note compressor models)
3. Oxus will return refurbished compressors with all accessories included and with like # of stroke (0.38" vs 0.44") as received:

0.44" / 9mm

0.44" / 10mm

0.38" / 9mm

0.38" / 10mm


Thomas 2450 Compressor (new)

Thomas 2450 Compressor (115V/60Hz)

Model # 2450AE44

RS-03020-1 (Individual Pack)

RS-03020-80 (80 Pack)



Thomas 2450 Compressor (230V/50Hz)

Model # 2450AUU44

RS-03027-1 (Individual Pack)

RS-03027-80 (80 Pack)

The new Thomas 2450 Compressors include:

0.44" Stroke Wobble Piston Compressor

Compressor Intake Tube

Capacitor

Great solution when you do not have a core; no core required.

Thomas 2450 Compressor Accessory Kits

RS-03030-1 (Compressor Intake Manifold Kit)

RS-03031-1 (Compressor Output Tube Kit)

RS-03032-1 (Compressor Spring/Washer Kit)

Oxus America

175 Kay Industrial Dr.

Orion Twp., MI 48359

Phone: 248-475-0925 | Toll Free: 888-475-1568

sales@oxusamerica.com

www.oxusamerica.com



Oxus is an engineering, manufacturing, and service company specializing in medical devices and gas separation technology. We are a FDA registered, MDSAP Certified, ISO 13485 Certified OEM Manufacturer and have been in business for over 20 years.