



**Oxus America**  
175 Kay Industrial Dr.  
Orion Twp., MI 48359  
Phone: 248-475-0925  
Toll Free: 888-475-1568  
sales@oxusamerica.com  
www.oxusamerica.com

## Instructions for the Use and Return of Stationary Sieve Cores - AirSep

### *Please read the following instructions prior to use:*

- *Each Oxus sieve is leak and performance tested to ensure purity at or above OEM specifications.*
- *Minimize the amount of time the sieve is exposed to the environment once the protective caps have been removed.*
- *Do not cut protective caps from sieves to avoid damaging fitting(s).*
- *Both sieve beds (“Right” and “Left”) must be replaced at the same time.*
- *Save the blue shipping cap to return with your core.*

### **Oxus sieves require return of an acceptable core within 3 weeks of shipment:**

- *Cores consist of any original OEM parts and/or any out-of-warranty Oxus refurbished parts only. Damaged or unusable cores will not be counted towards core receipt.*
- *Sieves with broken cap fittings will not be accepted as cores. Please do not return these sieves as cores.*
- *Return sieve cores in the same packaging as sent by Oxus:*
  - *Individual box with foam (2)*
  - *5 set pack box with inserts (2) and divider*
- *Sieves should be inserted in packaging for return shipment in the same manner as shipped from Oxus*
  - *Shipping cap facing down and the brass tee facing the center line of the box divider). “Right” assemblies to be placed in front row and “Left” assemblies to be placed in back row*
- *Ship core returns to: Oxus America, 175 Kay Industrial Dr. Orion Twp., MI 48359*



**For more information on our  
products and services:**

**Scan the QR Code  
OR**

**Visit us at [www.oxusamerica.com](http://www.oxusamerica.com)**



Oxus America  
175 Kay Industrial Dr.  
Orion Twp., MI 48359  
Phone: 248-475-0925  
Toll Free: 888-475-1568  
sales@oxusamerica.com  
www.oxusamerica.com

## Instructions for the Warranty Return of Stationary Sieve Cores - AirSep

*An RMA (Return Material Authorization) number is required prior to return of the product to Oxus for repair, credit, or replacement.*

- *All Oxus stationary sieves come with a 2 year warranty against leaks and rupture.*
- *Prior to return of the warranted product, contact Oxus through our designated RMA portal at [www.oxusamerica.com](http://www.oxusamerica.com) to initiate an RMA (must include serial number(s) and description of the problem).*
- *If applicable, Oxus will generate and provide an RMA number based on the ship date. The Oxus RMA number represents acknowledgment for return and tracking of product and does not guarantee warranty replacement.*
- *All RMA sieves **MUST** be capped immediately (both on the top and bottom fittings). **RMA sieves returned without caps will not be covered under warranty.***
- *RMA's should be shipped in a manner which provides protection. Sieves damaged in shipment due to improper packaging will not be covered under warranty.*
- ***RMA return sieves must be sent separately from return core shipments.***
- *The RMA number must be included on any shipping documentation and on the outside of the shipping container.*
- *Ship RMAs to: Oxus America, 175 Kay Industrial Dr. Orion Twp., MI 48359*

---

### Also available from Oxus...

*Stationary Sieve Replacements for Invacare, Respironics and DeVilbiss*

*POC Sieve Replacements for AirSep Freestyle 3 & 5*

*Stationary Valve Replacements for Invacare and Respironics EverFlo*

*Compressor Rebuild Kits*