



SERVICE SOLUTIONS
Stationary Concentrator
Valves



Oxus has developed a high quality aftermarket solution for Stationary Concentrator replacement valves. Oxus valves are designed to meet or exceed OEM specifications and come with a 1 year warranty against wear and workmanship.

Oxus Refurbished Valves for Invacare and DeVilbiss Concentrators:

All Oxus refurbished valves include:

- Complete disassembly and ultrasonic cleaning of all valve components
- DeVilbiss Rotary Valve refurbishment also includes 100% replacement of top caps and sealing/sliding interfaces, with new parts
- Lubrication using air-path compatible lubricant
- 100% functionally tested to confirm valve performance prior to shipment

Oxus Valves for Respironics EverFlo Concentrators:

- 100 % new products
- Designed to eliminate body cracking

Oxus America
175 Kay Industrial Dr.
Orion Twp., MI 48359
Phone: 248-475-0925 | Toll Free: 888-475-1568
sales@oxusamerica.com
www.oxusamerica.com



Oxus is an engineering, manufacturing, and service company specializing in medical devices and gas separation technology. We are a FDA registered, MDSAP Certified, ISO 13485 Certified OEM Manufacturer and have been in business for over 20 years.

APPLICATION	OEM P/N	DESCRIPTION	OXUS P/N
All Invacare Concentrators	1101141	 4-Way Valve	RS-04000-1 (individual)
			RS-04000-24 (24 Pack)
Respironics EverFlo (Pre-June 2013)	1111248	 Valve Body - Internal Pilot - SMC Replacement	RS-04011-1 (individual)
			RS-04011-25 (25 Pack)
Respironics EverFlo (Pre-June 2013)	1114952	 Valve Body - Internal Pilot - ASCO Replacement	RS-04012-1 (individual)
			RS-04012-25 (25 Pack)
Respironics EverFlo (Post-June 2013)	1115003	 Valve Body - External Pilot - ASCO Replacement	RS-04013-1 (individual)
			RS-04013-25 (25 Pack)
All Respironics EverFlo Concentrators	1111249	 SMC Pilot Valve Asm.	RS-04004-1 (individual)
			RS-04004-25 (25 Pack)
All Respironics EverFlo (ASCO Applications)	1114954	 Pilot Valve Asm. - ASCO Replacement	RS-04005-1 (individual)
			RS-04005-25 (25 Pack)
Drive/DeVilbiss 5L (525DS)	515ADZ-702 515D-702	 5L Rotary Valve	RS-04050-1 (individual)
			RS-04050-20 (20 Pack)
Drive/DeVilbiss 10L (1025DS)	1025D-702	 10L Rotary Valve	RS-04051-1 (individual)
			RS-04051-20 (20 Pack)

Oxus America
175 Kay Industrial Dr.
Orion Twp., MI 48359
Phone: 248-475-0925 | Toll Free: 888-475-1568
sales@oxusamerica.com
www.oxusamerica.com



Oxus is an engineering, manufacturing, and service company specializing in medical devices and gas separation technology. We are a FDA registered, MDSAP Certified, ISO 13485 Certified OEM Manufacturer and have been in business for over 20 years.