



SERVICE SOLUTIONS

Stationary Sieve Replacements




Oxus has developed a high quality aftermarket solution for Stationary Concentrator replacement sieves with an eye toward cost reduction.

All Oxus Stationary Concentrator replacement sieves:

- 100% leak and performance tested to ensure purity at or above OEM standards, ensuring full compliance with FDA guidelines
- Critical parts mechanically cleaned
- Maximum zeolite fill to ensure long-term durability
- **2 year warranty against leaks and rupture**

APPLICATION	OEM PART #	DESCRIPTION	OXUS PART #
 Invacare Perfecto2V	1165099	 Perfecto2V Sieve	RS-02002
 Invacare Perfecto2	1161257	 Perfecto2 Sieve	RS-02003
 Invacare Platinum 9/10	1144034	 Platinum 9/10 Sieve	RS-02000
 Invacare Platinum 5	1101707 1161699	 Platinum 5 Sieve	RS-02001

APPLICATION	OEM PART #	DESCRIPTION	OXUS PART #
 Respironics EverFlo	1038825 0.035 x 11.5"	 EverFlo Sieve (0.035 x 11.5")	RS-02011 RS-02014 (export)
	1052382 0.032 x 12.5"	 EverFlo Sieve (0.032 x 12.5")	RS-02012 RS-02015 (export)
	1050248 0.035 x 12.5"	 EverFlo Sieve (0.035 x 12.5")	RS-02013 RS-02016 (export)
 DeVilbiss 5L (525DS)	525D-619	 DeVilbiss 5L Sieve standard and long versions available	RS-02020 (Standard) RS-02021 (Long)
 DeVilbiss 10L (1025DS)	1025D-619	 DeVilbiss 10L Sieve standard and long versions available	RS-02030 (Standard) RS-02031 (Long)
 Caire Newlife Elite 5L	BE001-7R	 Caire AirSep Newlife Elite Sieve	RS-02065
 Respironics Millennium M10	1017774	 Millennium M10 Sieve	RS-02070
 Direct Supply Attendant 5L	FYM38	 Direct Supply Attendant 5L Sieve	RS-02080
 React Stratus 5	OEM P/N Not Available	 Stratus 5L Sieve	RS-02090